



# The specific case

Surgical treatment  
of a olecranon fracture,  
partial non-union with  
**CERASORB® Foam**

Drs. E.J. Hekma  
Rijnstate Hospital  
Arnhem, NL

**CERASORB®-PROMISE**

**CERASORB®.**

Keeps its words in bone regeneration.

## The specific case

# + Olecranon fracture, partial non-union

## Surgical treatment of a olecranon fracture, partial non-union with CERASORB® Foam

*„In this case, the use of the synthetic, fully resorbable pure-phase  $\beta$  - Tricalcium-phosphate **CERASORB®** is an appropriate solution for the management of a partial non-union. Moreover, the foam version of **CERASORB®** has shown excellent intraoperative handling characteristics with an unparalleled ease-of-use.“*

Drs. E.J. Hekma - Trauma Surgeon  
Rijnstate Hospital  
Arnhem, NL

[Fig. 1]



[Fig. 2]



### Anamnesis

A 57-year old male crashed with his bicycle resulting in a comminuted right olecranon fracture [Fig. 1]. Two superficial lesions could be observed on his elbow. He was operated on November 12, 2014 by using internal plate fixation via ORIF (Synthes® LCP Olecranon plate)

T = Months

3

In February 2015 the patient showed signs of an infected osteosynthesis. Osteolysis around one of the screws was suspected. The patient was treated with Augmentin®.


No signs of consolidation could be observed in the follow-up period [Fig. 3] and the patient was re-operated on June 26, 2015.

[Fig. 3]



A new plate was implanted following debridement and the bony defect was grafted with autologous bone chips harvested from the tibial head [Fig. 4]. The treatment with Augmentin® was continued.

[Fig. 4]




In October 2015, the patient underwent electrostimulation to stimulate bone healing, however, without any satisfying result.

### Surgical treatment with CERASORB® Foam

On January 25, 2016 the patient was operated for the third time. X-ray and CT scans had revealed a relatively large partial non-union [Figs. 5 and 6]. There were no signs of infection.

[Fig. 5]



[Fig. 6]



The plate was removed and the defect was debrided and filled with 10 cc of **CERASORB® Foam** (synthetic, fully resorbable, pure-phase  $\beta$ -TCP + porcine collagen matrix bone graft substitute) [Fig. 7].

The patient received a removable cast for 6 weeks. Directly post-op he started with isometric exercises.

[Fig. 7]



#### Follow-up and discharge

The x-ray of March 4, 2016 was indicative of an ongoing remodeling process, the bone graft substitute being gradually resorbed and replaced by endogenous new bone [Fig. 8].

[Fig. 8]



On June 6, 2016, the x-ray showed a far progressed remodeling into new bone [Fig. 9].

Any **CERASORB® Foam** present in small amount in the soft tissues was resorbed, with no risk of ectopic bone formation. The patient had a full range of motion, no more complaints, was back to work and had restarted cycling.

[Fig. 9]



16

19

## Take home messages

- + Grafting with **CERASORB® Foam** has shown to be an appropriate solution in this case.
- + **CERASORB® Foam**, a fully resorbable, synthetic pure-phase  $\beta$ -TCP-collagen bone graft substitute, can offer appropriate solutions for the management of partial non-unions.
- + It can be discussed whether the addition of bone marrow aspirate or autologous bone graft material to the osteoconductive scaffold **CERASORB® Foam** may be necessary to manage these kinds of defects which are at higher risk of non-healing (partial non-union, specific drug usage, low grade infection in anamnesis). In this case, bone marrow aspirate was not added, still resulting in the desired outcome.

## + Maximum flexibility

### CERASORB® Foam



**CERASORB® Foam** is a pure-phase, fully synthetic  $\beta$ -TCP-collagen scaffold. The product is available in both, a mouldable and flexible foam version.

### Type of use:

- + As stand alone in critical-size metaphyseal or segmental defects not requiring additional cells and growth factors.
- + In combination with autologous materials in defects requiring additional cells and growth factors.
- + In combination with autologous bone chips
- + In combination with RIA
- + In combination with bone marrow aspirate

### Handling:

- + **CERASORB® Foam** can be moulded when soaked with patients blood or bone marrow aspirate in a ratio of 1:1
- + **CERASORB® Foam** can be cut in small pieces to facilitate combination with autologous materials.

**CERASORB®** bone-regeneration materials.  
We offer tailor-made solutions for diverse requirements.

**You have our word!**

**curasan**  
Regenerative Medicine

curasan AG  
Lindigstrasse 4  
63801 Kleinostheim  
Germany  
cerasorb@curasan.de  
www.curasan.de



**CERASORB® – Always First Choice.**